

MICRO SURGERY

[MH18-867](#) BISHOP-HARMON Micro Tissue Forceps, 3-3/8" (8.6 cm), delicate 0.3 mm 1 X 2 teeth

[MH18-835](#) BONN Micro Suturing Forceps, 2-3/4" (7 cm), 1 X 2 teeth, 0.12 mm, with tying platform

[MH18-836](#) BONN Micro Suturing Forceps, 3-3/4" (9.5 cm), 1 X 2 teeth, 0.12 mm, with tying platform

MH17-2602 JACOBSON Micro Mosquito Forceps, curved, extremely delicate

[MH17-2600](#) JACOBSON Micro Mosquito Forceps, straight, extremely delicate

[MH18-949](#) MCPHERSON Micro Iris Suturing Forceps, 3-1/2" (8.9 cm), angled, smooth platform 0.3 X 5 mm, no teeth

[MH18-948](#) MCPHERSON Micro Iris Suturing Forceps, 3-1/2" (8.9 cm), straight, smooth platform 0.3 mm X 5 mm, no teeth

MH17-2510 Micro ADSON Dressing Forceps, 4-3/4" (12.1 cm), serrated jaws, 0.5 mm wide

[MH17-2500](#) Micro ADSON Tissue Forceps, 4-3/4" (12.1 cm), tips 0.5 mm wide, 1 X 2 teeth

[MH19-2110](#) Micro Ear Forceps, Alligator Type, 3-1/4" (8.3 cm) shaft, oval cup jaws, 0.8 mm wide

[MH19-2081](#) Micro Ear Forceps, Alligator Type, 3-1/4" (8.3 cm) shaft, serrated jaws, 0.6 mm wide