

NEUROSURGICAL

- [MH6-198](#) ADSON Bayonet Dressing Forceps 8-1/4" (21 cm), serrated
- MH6-118 ADSON Dressing Forceps, 4-3/4" (12.1 cm), delicate, serrated
- [MH6-123](#) ADSON Dressing, Tissue " Suture Forceps 1 X 2 teeth, 4-3/4" (12.1 cm), with tying platform
- [MH26-1090](#) ADSON Dura and Dissecting Hook, 8" (20.3 cm), sharp blade 5 mm at right angle
- [MH7-242](#) ADSON Forceps, 7-1/4" (18.4 cm), curved
- [MH7-240](#) ADSON Forceps, 7-1/4" (18.4 cm), straight
- [MH11-702](#) ADSON Retractor, 7-1/2" (19.1 cm), sharp, 4 X 4 teeth 3/4" (1.9 cm) X 3/4" (1.9 cm), angled arms
- MH11-700 ADSON Retractor, 7-1/2" (19.1 cm), sharp, 4 X 4 teeth 3/4" (1.9 cm) X 3/4" (1.9 cm), straight Arms
- [MH6-121](#) ADSON Tissue Forceps 1 X 2 teeth, 4-3/4" (12.1 cm), cross serrated tips, delicate
- [MH6-120](#) ADSON Tissue Forceps 1 X 2 teeth, 4-3/4" (12.1 cm), delicate
- [MH6-122](#) ADSON Tissue Forceps 2 X 3 teeth, 4-3/4" (12.1 cm), delicate
- [MH19-580](#) BARON Suction Tube, 3 french (1 mm), with finger cut-off, working length 7.5 cm
- [MH19-582](#) BARON Suction Tube, 5 french (1.7 mm), with finger cut-off, working length 7.5 cm
- [MH19-584](#) BARON Suction Tube, 7 french (2.3 mm), with finger cut-off, working length 7.5 cm
- [MH11-684](#) BECKMAN-ADSON Retractor, 12-1/2" (31.8 cm), hinged, sharp, 4 X 4 teeth, 3/4" (1.9 cm) X 1" (2.5 cm)
- MH11-630-BL BECKMAN-WEITLANER Retractor, 5-1/2" (14 cm), blunt, 3 X 4 teeth, with hinge blades
- MH11-630-SH BECKMAN-WEITLANER Retractor, 5-1/2" (14 cm), sharp, 3 X 4 teeth, with hinge blades
- [MH19-844](#) BEYER Rongeur, 7" (17.8 cm), double action, slightly curved jaws 3.5 mm wide
- [MH6-124](#) BROWN-ADSON Tissue Forceps, 7 X 7 teeth, 4-3/4" (12.1 cm)
- [MH6-192](#) CUSHING Dressing Forceps 1 X 2 teeth, 7-1/4" (18.4 cm), bayonet, scraper end
- [MH6-190](#) CUSHING Dressing Forceps, 7-1/4" (18.4 cm), bayonet, serrated, scraper end
- [MH26-394](#) CUSHING Pituitary Rongeur, 2 X 10 mm cup jaws, 5" (12.7 cm) shank, straight
- MH26-397 CUSHING Pituitary Rongeur, 2 X 10 mm cup jaws, 7" (17.8 cm) shank, straight
- [MH11-102](#) CUSHING Vein Retractor, 9" (22.9 cm), blade 11 X 13 mm
- [MH7-102](#) DANDY Forceps, 5-1/2" (14 cm), curved sideways
- [MH26-956](#) FERRIS-SMITH Tissue Forceps, 7" (17.8 cm), serrated platform tips with 1 X 2 teeth
- MH26-958 FERRIS-SMITH Tissue Forceps, 7" (17.8 cm), serrated platform tips with 2 X 3 teeth
- [MH26-1084](#) FRAZIER Dura Hook, sharp point
- MH19-575 FRAZIER Suction Tube, 10 french, 3.3 mm, with finger cut-off, angled
- MH19-576 FRAZIER Suction Tube, 11 french, 3.6 mm, with finger cut-off, angled
- MH19-577 FRAZIER Suction Tube, 12 french, 4 mm, with finger cut-off, angled
- MH19-570 FRAZIER Suction Tube, 6 french, 2 mm, with finger cut-off, angled
- MH19-572 FRAZIER Suction Tube, 7 french, 2.3 mm, with finger cut-off, angled
- MH19-573 FRAZIER Suction Tube, 8 french, 2.6 mm, with finger cut-off, angled
- MH19-574 FRAZIER Suction Tube, 9 french, 3 mm, with finger cut-off, angled
- [MH18-1968](#) FREER Elevator, 7" (17.8 cm), double end sharp and blunt blades 5 mm wide
- MH11-364 GELPI Retractor 3.5 (8.9cm), with grip lock, Mini
- MH11-366 GELPI Retractor, 7-1/2" (19.1 cm), with grip lock
- [MH11-365](#) GELPI Retractor, Pediatric, 5-1/2" (14 cm), with grip lock
- [MH6-180](#) GERALD Dressing Forceps, 7" (17.8 cm), straight, serrated
- [MH6-182](#) GERALD Tissue Forceps 1 X 2 teeth, 7" (17.8 cm), straight
- [MH26-124](#) GIGLI-STRULLY Saw Handle, standard pattern, loop grip
- [MH26-119](#) GIGLI Saw Handle, standard pattern, T-shaped, turnable
- [MH6-196](#) GRUENWALD (ADSON) Forceps, 8.25" (21 cm)
- [MH6-197](#) GRUENWALD (JANSEN) Bayonet Dressing Forceps, 6-1/2" (16.5 cm), serrated