

SCALPEL HANDLES

[MH4-10](#) Knife Handle no. 3L, for deep surgery, fitting surgical blades nos. 10 thru 15c, extra fine.

[MH4-7](#) Knife Handle no. 3, fitting surgical blades nos. 10 thru 15c, extra fine

[MH4-8](#) Knife Handle no. 4, fitting surgical blades nos. 20 thru 25, extra fine

[MH4-9](#) Knife Handle no. 7, fitting surgical blades nos. 10 thru 15c, extra fine